


Specifications Sheet	No.SPEC-3306E	(1st Edition)	
<b>OM7000A</b>	<b>COULOMETER</b>		

## OUTLINE

The OM7000 is a coulometric bod measuring equipment developed by government industrial department laboratory, Hokkaido and Ohkura. Its application extends to measuring decomposition of chemicals by bacteria (microscopic organism), sewage disposal simulation in waste treatment, automating and quantifying tissue culture in pathology. Evaluating effectiveness of vaccine and antibiotics, measurement of respiratory and oxygen consumption in chemical, marine products and agricultural industries.

Can perform the Print of various calculations and the test result in combination with data processing software by having changed the Measuring unit from the paper chart Expression recorder to the touch panel system indicator.

## FEATURES

- Measurement by automated oxygen supply  
Be most suitable for straight Degradating test such as OECD test guidelines 301C, JIS K 6950.
- Touch panel system indicator  
Adopted an exclusive touch panel system indicator as an operation part and display. Can confirm a result as a chart by simple operation and a trendy display.
- Be possible in combination with conventional products  
The Measuring unit is accessible to conventional thermostatic ovens such as OM7000. Can just use the glasswares such as the Electrolytic bottle, too.
- Data processing software  
Prepared data processing software for Windows7. Can perform processing such as a calculation, the Print of data by being able to be connected to eight Measuring unit with one PC, and using data processing software.
- Estimates BOD<sub>5</sub> using raw sample solution
- One time measurement of 6 samples in temperature bath
- Direct readout of oxygen consumption in "ppm" or "mg"
- O-ring sealed cultivating bottle




## BLOCK DIAGRAM

The OM7000A consists of a detector, a temperature bath and a measuring unit.

This Coulometer is comprised of a detecting element, a thermostatic oven part, Measuring unit. In addition, the Coulometer data processing section is comprised of a PC, data processing software (for exclusive use of the Coulometer) and an LAN cable.

The private telecommunication boards are still less necessary because they perform the communication between Measuring unit - PC in an Ethernet.



The PC body does not include in annotation.

## SPECIFICATIONS

**Sensitivity:** 0.2 ppm

**Operation temperature-humidity:** 5~30°C, 35 - 85%RH

**Power input Voltage:** 100VAC±10VAC 50/60Hz,  
220VAC±22VAC 50/60Hz

**Power consumption:** Approx. 1070VA

**Equipment weight:** Approx. 220kg

## Measuring unit

**Measurement range:** 15 mg (50ppm), 30 mg (100ppm),  
75 mg (250ppm), 150 mg (500ppm), five phases of  
reshuffling of 300 mg (1,000ppm)

**Measurement score:** 6 points

**Temperature input:** Score, 2 points

**Temperature element,;** Pt100

**Measuring range,;** 0.0-100.0°C

**Accuracy: Electrolytic current;** ±1% of rdg

**Temperature input;** ±0.5°C

**Indicator:** 5.7 inches of TFT color LCD (320\*240 dot)  
touch panels, with backlight. Display measured value,  
a chart.

**Operation button: The number of the buttons ;** 3

**Function,** REC, Test initiation / stop

MENU, Various setting windows display

FUNC, Do not use

**Record function: Outside recording medium;**

SD memory card (available only for an operation  
check product correspondence in SD/SDHC technical  
standard)

internal memory, Approximately 100MB file  
preservation period, One minute

**Alarm function:** Warning kind, Temperature warning

•High-low-limit check alarm

•Deviation alarm

Totalizing monitoring warning

•Manometer continuation ON

•Manometer continuation OFF

•Integral velocity abnormality

•Electrolysis current value abnormality

**Alarm output:** Output, 2 points

Contact rating, 3A/250VAC, 3A/30VDC

**Communication:** An Ethernet, a transmission type:  
10BASE-T

**Memory back up:** Test data and the parameter are the  
inside. Save in flash memory. The clock backs up with  
a built-in lithium battery (battery life approximately five  
years at the time of no conducting).

**Power consumption:** Approx. 55VA Max

**Weight:** Approx. 5 kg

## DETECTOR

**Manometer:** Sensitivity 9.806Pa (1mAq, 0.2 ppm), silver  
Electrode, contents of copper sulfate  
solution

**Electrolytic bottle:** Capacity 500ml, electrode platinum  
and copper, contents saturated copper sulfate solution  
and 5% sulfuric acid

**Cultivating bottle:** 300ml with absorbent container

**Stirrer:** 540rpm/50Hz or 650rpm/60Hz for 6 bottles

## TEMPERATURE BATH

**Method:** Forced ventilation Expression

**Dimension:** W1050 X H1363 X D710mm

**Ambient air temperature:** 5~30°C

**Control mode:**

SSR drive, ON/OFF control (refrigerator continuous  
running)

**Temperature set point:** 5~40°C

**Temperature distribution:** ±1°C

**Power supply:** 100VAC 50/60Hz or 220VAC 50/60Hz  
(the deployment trans other than 100VAC)

**Power consumption:** Approx. 1kVA

**Equipment weight:** Approx. 180 kg

**Refrigerator:** Reciprocating type refrigerator (180W)  
Refrigerant, R-134a

**Temperature controller:**

Digital temperature controller made in Ohkura Electric  
Type, EC5502R00020

-150.0~+150.0°C, Pt100

**Magnetic Stirrer specifications:**

Permanent magneto Magnetic Stirrer, Six

(1) Standard specifications

100VAC 3W 540rpm/650rpm turns clockwise  
(50/60Hz)

(2) Direction of rotation inversion specifications

With timer for the Magnetic Stirrer inversion (all  
channels are common specifications and  
equivalence standard as for the revolution  
speed) Turn clockwise, and the timer 1~  
60sec setting for the timer 1 ~10min setting  
possible stop for rotating anticlockwise is  
possible

(3) Revolution speed variableness specifications

Setting possibility range 70~800rpm available  
for revolution speed setting every each channel

**Thermostatic oven power cable:**

With 100VAC specifications 3P plug

## Data processing system

**Operating environment:** OS, Windows7

**Junction number:** Up to eight (for none of the PCs)

**Junction method:** Ethernet

## STANDARD ACCESSARIES



DESCRIPTION	QUANTIFY	DESCRIPTION	QUANTIFY
Manometers including spares	7	Silicone grease (50g)	1
Cultivating bottle including spares	7	Seal tape	1
Air tank including spares	7	Liquid paraffin	1
Electrolytic bottle kit	6	Beaker (1l for drain)	1
Tube	18	Fuse (15A) for heater	1
Cable	1	SD memory card (SLC type)	1
Manometer rubber cap with electrode	20	Coulometer instruction manual	1
Stirrer chip (Φ10 ×50)	6	Temperature controller instruction manual	1
Copper sulfate (500g)	2		
Soda lime (500g)	1		

## MODEL CODE NUMBER



MODEL	POWER INPUT	STIRRER	Configu ration	Option	Software	Special specifications	DESCRIPTION
OM7000A							
	1						100VAC 50/60 Hz
	2						220VAC 50/60 Hz
	9						Custom (consult factory)
		0					None
		1					Standard single direction cw, 540 rpm at 50Hz, 650rpm at 60Hz
		2					Bi – directional, 10 min. by cw, 1 min. stop and 10 min. by ccw, 540rpm/ 50Hz 650/ 60Hz
		3					Single direction cw, adjustable rotation between 70 and 800 rpm
			1				Complete set
			2				Only as for the Measuring unit
				0			None
				9			Custom (consult factory)
					0		None
					1		Standard specifications
					2		Straight Degradating test
					9		Others
						0	None
						9	Custom (consult factory)

## DIMENSIONS (Unit: mm)

### ●MEASURING UNIT



### ●TEMPERATURE BATH



Note 1) Windows7, Excel are the registered trademarks of Microsoft Corporation of the U.S.A.  
Note 2) Ethernet is the registered trademark of Xerox Corporation.

## ⚠ CAUTION

Do not install this device before consulting instruction manual

**Ohkura**

OHKURA ELECTRIC CO., LTD.

Head Office / Factory

Sales Offices

URL

e-mail (in English)

Saitama, JAPAN

Tokyo, Osaka, Nagoya, Kyushu

<http://www.ohkura.co.jp/>

[contact\\_e@ohkura.co.jp](mailto:contact_e@ohkura.co.jp)

Specifications are subject to change without notice.

For further information, a quotation or a demonstration please contact to:

Printed in Japan: June 2016